

24.0 EXHIBIT 15.0 ESD GLOVE FORM (EXAMPLE)

ANSI/ESD S20.20-2014 Compliance Verification													
ESD Glove													
Area:													
Date:													
ESD Auditor:													
		Temp (°F/°C):											
		Humidity (RH%):											
	Surface Resistance Test	Surface Resistance Test	Surface Resistance Test	Surface Resistance Test	Surface Resistance Test	Surface Resistance Test	Comments/Issues (Actions Taken)						
Finger		Sample 1	Sample 2	Sample 3	Sample 4	Sample 5							
	(No Glove)	(With Glove)	(With Glove)	(With Glove)	(With Glove)	(With Glove)							
Rg < 1 x 10 ⁹ Ohms													
	Actual Measurement	P/F	Actual Measurement	P/F	Actual Measurement	P/F	Actual Measurement	P/F	Actual Measurement	P/F	Actual Measurement	P/F	
1 st finger index													
2 nd finger													
3 rd finger													
4 th finger													
5 th finger thumb													
Results		Test Equip: PRS-801 Prostat Meter											
<10 ⁸		White											
>10 ⁸ <10 ⁹ Process Indicator		Yellow											
>10 ⁹ Fail		Red											
		PCF-825 Prostat Café Test Fixture											
		PWS-611-PGC Prostat No Resistor Ground Cord (RED)											
		Alcohol Wipe - Electrode Cleaning Wipe											
		DIST: QAS											
QA178Log-Glove Rev C													

25.0 EXHIBIT 16.0 ESD BAG & FOAM FORM (EXAMPLE)

ANSI/ESD S20.20-2014 Compliance Verification												
ESD Bag & Foam												
Area:												
Date:												
ESD Auditor:												
		Temp (°F/°C):										
		Humidity (RH%):										
Item ID	Surface Resistance Test	Surface Resistance Test	Comments/Issues (Actions Taken)									
	Rg < 1x10 ¹¹ Ohms	Rg < 1x10 ¹² Ohms per sq inch										
	(Increase Exponent by 1 for ohm/sq inch calc.)											
	Actual Measurement	P/F	Actual Measurement	P/F								
Results		Test Equip: PRS-801 Prostat Meter										
<10 ¹¹		White										
>10 ¹¹ <10 ¹² Process Indicator		Yellow										
>10 ¹² Fail		Red										
		PRF-911 Prostat Concentric Ring Test Fixture										
		PTB-920 Prostat Dual Function Test Bed										
		Cleaning Wipe - Electrode Cleaning Wipe										
		DIST: QAS										
QA178Log-BagFoam Rev C												